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0.1 → 50

Macro zoom lens

With a range from 0.1x - 50x magnification, a target can be viewed from its entirety down to more in-depth observation. This macro lens excels in workability and high performance with click-style magnification adjustment, an aperture mechanism, and a viewing distance of 95 mm or more.

Model		VH-Z00R/Z00T						
Magnification		0.1x	0.5x	1x	5x	10x	30x	50x
Field-of-view (mm)	Horizontal	3200	640	320	61	30.5	10.2	6.1
	Vertical	2400	480	240	45.5	22.8	7.6	4.6
	Diagonal	4000	800	400	76.2	38.1	12.7	7.6
Working distance (mm)		Approx 7700	Approx 1500	Approx 720	95			

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20 → 200

Versatile lens provides high-resolution imaging with large depth-of-field

The VH-Z20R/Z20T offers high-resolution observation at general purpose magnifications of 20x - 200x. This lens has been designed to optimise both depth-of-field and resolution and can also be used in handheld mode.

Model		VH-Z20R/Z20T						
Magnification		20x	30x	50x	100x	150x	200x	
Field-of-view (mm)	Horizontal	15.24	10.16	6.10	3.05	2.03	1.52	
	Vertical	11.40	7.60	4.56	2.28	1.52	1.14	
	Diagonal	19.05	12.70	7.62	3.81	2.54	1.91	
Depth-of-field ² (mm)		34	15.5	6.0	1.6	0.74	0.44	
Working distance (mm)		25.5						

1. Magnification on a 15-inch monitor

2. The value when the lens is set with priority to depth-of-field. The depth-of-field changes depending on the setting of the aperture.

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VH-Z100R/Z100T

100 → 1000

High-performance lens with long working distance

This innovative lens was developed to satisfy the need for high-resolution, long working distance, and large depth-of-field. Provides both ring light and bright field illumination.

Model		VH-Z100R/Z100T						
Magnification		100x	200x	300x	500x	700x	1000x	
Field-of-view (mm)	Horizontal	3.05	1.53	1.02	0.61	0.44	0.30	
	Vertical	2.28	1.14	0.76	0.46	0.33	0.23	
	Diagonal	3.81	1.90	1.27	0.76	0.54	0.38	
Working distance (mm)		25 (20 ²)						

1. Magnification on a 15-inch monitor

2. When the triple illumination adapter is attached.

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NEW

250 → 2500

Observe with both bright field and dark field at high-magnification

Easily switch between ring light and coaxial illumination with just the touch of a button. View objects at up to 2500x magnification while still maintaining a 6.5 mm working distance.

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500 → 5000

Our highest magnification/resolution zoom lens

This zoom lens incorporates high-quality fluorite optics to provide the highest resolution in its class. With an N.A. of 0.82, achieve up to 5000x magnification with a 4.4 mm working distance.

Model		VH-Z500R/Z500T				
Magnification		500x	1000x	2000x	3000x	5000x
Field-of-view (mm)	Horizontal	610	305	152	102	61
	Vertical	457	229	114	76	46
	Diagonal	762	381	191	127	76
Working distance (mm)		4.4				

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VH-Z20UR/Z20UT **NEW**

20 200

Optimal lighting with the touch of a button

This newly-designed lens has the ability to perform bright/dark field and DIC observation, even at lower magnification ranges. A unique illumination system allows users to switch between three different types of lighting by simply pressing a button.

Bright-field	Dark-field
Partial	DIC

Model		VH-Z20UR/Z20UT					
Magnification		20x	40x	80x	100x	160x	200x
Field-of-view (mm)	Horizontal	15.24	7.62	3.81	3.05	1.91	1.52
	Vertical	11.40	5.70	2.85	2.28	1.43	1.14
	Diagonal	19.05	9.53	4.76	3.81	2.38	1.91
Working distance (mm)		20.8 ²					

1. Magnification on a 15-inch monitor.
2. With the wide-area illumination attachment equipped.

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VH-Z100UR/Z100UT

100 1000

Differential Interference Contrast (DIC) lens


Bright/dark field, polarised transmitted, and DIC observation can be performed with this lens. DIC observation makes it possible to clearly visualise surface topography of low-contrast and transparent objects - typically difficult with conventional bright field lighting.

Bright-field	Dark-field
Polarisation	DIC

Model		VH-Z100UR/Z100UT					
Magnification		100x	200x	300x	500x	700x	1000x
Field-of-view (mm)	Horizontal	3.05	1.53	1.02	0.61	0.44	0.30
	Vertical	2.28	1.14	0.76	0.46	0.33	0.23
	Diagonal	3.81	1.90	1.27	0.76	0.54	0.38
Working distance (mm)		25(20 ²)					

1. Magnification on a 15-inch monitor.
2. When the triple illumination adapter is attached.

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50 500


Long Range Lens with a 85 mm Working Distance

Enables high-magnification observation while maintaining a long working distance. This lens is ideal for viewing objects that have highly-irregular surfaces or recesses that cannot be observed up close.

Model		VH-Z50L/Z50T					
Magnification		50x	100x	200x	300x	400x	500x
Field-of-view (µm)	Horizontal	6.09	3.05	1.53	1.02	0.76	0.61
	Vertical	4.57	2.28	1.14	0.76	0.57	0.46
	Diagonal	7.62	3.81	1.90	1.27	0.95	0.76
Working distance (mm)		85					

1. Magnification on a 15-inch monitor.

VHX Series System Line Up

Lens bracket
**VHX-J00T/
 J20T/J20UT/
 J50T/J100T/
 J250T/J500T**

Bayonet-type
 standard attachment
**OP-51478/
 OP-66871** ⁶

VHX-5100

VHX-5020

VHX-5000

VHX-1020

VHX-700FE


Bayonet-type C-mount
 attachment
OP-51479

Camera lens (C-mount type)


Microscope

VH-F

VHX-H2M2 ⁸
 XY measurement
 software


VHX-H4M ⁹
 3D measurement
 software

PC

PC

PC

PC

Communication software
 (2D measurement,
 3D/HDR display software)

LAN

CD-R/RW

USB memory, etc.

USB

Hard disk drive

DVI

Monitor

A large screen monitor can be
 connected when several people need
 to view a target simultaneously.
 Since the target image is magnified
 further on the large screen, it can be
 checked more thoroughly.

3D Profile measurement unit
VHX-S15

Foot switch

OPTION

Free-angle
 observation system
 (Z-axis motorised,
 manual XY)
VH-X-S500E

Free-angle observation
 system (XYZ-axis motorised
 with transmitted illumination)
VH-X-S550E


High-precision stand
 (with X-Y stage and
 transmitted illumination)
VH-S5

VH lens mounting stand
 (with XY stage)
**OP-25539
 OP-22124**

X-Y measurement system
VH-M100E

Transmitted
 illumination unit
OP-84484

Display unit
OP-84483

1. TRIPLE-R compliant lenses VH-Z00T/Z20T/Z20UT/Z50T/Z100T/Z250T/Z500T are fitted with Automatic Lens/Zoom Recognition units and connection recognition mount, respectively.
2. OP-72407 and OP-72406 are required when coaxial illumination is used.
3. Included with the VH-Z20UR/Z20UT.
4. The optional bore fibre cable (OP-87201) is required.
5. The optional light guide attachment (either OP-51482 or OP-87790) is required.
6. OP-66871 is required when the VH-Z00R, Z20R, or Z25 is used.
7. A C-mount adapter suitable for the microscope is required.
8. VHX-H1M1 is required for the VHX-700FE.
9. VHX-H3M is required for the VHX-700FE.

SPECIFICATIONS

Basic functions: Controller

Model		VHX-5000	VHX-700FE	
Camera	Image receiving element	1/1.8-inch, CMOS image sensor Virtual pixels: 1600 (H) x 1200 (V)	1/1.8-inch, CCD image sensor Virtual pixels: 1600 (H) x 1200 (V)	
	Scan method	Progressive	Progressive	
	Frame rate	50 frames/sec. (max.)	15 frames/sec. and 28 frames/sec. selectable	
	Resolution	Normal	1600 (H) x 1200 (V) Approx. 1000 TV lines	1600 (H) x 1200 (V) Approx. 1000 TV lines
		3CMOS ^{1,3}	1600 (H) x 1200 (V) Approx. 1200 TV lines (2 million pixels x 3CMOS mode, Excellent colour reproducibility)	–
		High resolution ³	3200 (H) x 2400 (V) Approx. 1600 TV lines	–
		Super high resolution ³	4800 (H) x 3600 (V) Approx. 2000 TV lines or more	–
	Super high resolution x 3CMOS ^{2,3}	4800 (H) x 3600 (V) Approx. 2000 TV lines or more (18 million pixels x 3CMOS mode, Excellent colour reproducibility)	–	
	High Dynamic Range	16-bit resolution through RGB data from each pixel	–	
	Gain	AUTO, MANUAL, PRESET	AUTO, MANUAL, PRESET	
	Electronic shutter	AUTO, MANUAL, 1/60, 1/120, 1/250, 1/500, 1/1000, 1/2000, 1/5000, 1/9000, 1/19000	AUTO, MANUAL, 1/15, 1/30, 1/60, 1/120, 1/250, 1/500, 1/1000, 1/2000, 1/5000, 1/9000, 1/19000	
	Supercharge shutter	0.2 sec. to 4 sec. Can be set in increments of 0.1 sec.	0.2 sec. to 17 sec. Can be set in increments of 0.1 sec.	
	White balance	AUTO, MANUAL, ONE-PUSH SET, PRESET (2700K, 3200K, 5600K, 9000K)	AUTO, MANUAL, ONE-PUSH SET, PRESET (2700K, 3200K, 5600K, 9000K)	
Back-focus adjustment	Not required	Not required		
LCD monitor ⁵	Size	Colour LCD (IPS) 23"	Colour LCD (TFT) 17"	
	Panel size	509.184 (H) x 286.416 (V) mm	365.76 (H) x 228.60 (V) mm	
	Pixel pitch	0.2652 mm (H) x 0.2652 mm (V)	0.1905 mm (H) x 0.1905 mm (V)	
	Number of pixels	1920 (H) x 1080 (V) (FHD)	1920 (H) x 1200 (V) (WUXGA)	
	Display colour	Approx. 16,770,000 colours ⁴	Approx. 16,770,000 colours ⁴	
	Brightness	300 cd/m ² (Centre 1 Point, typical)	270 cd/m ² (typical)	
	Contrast ratio	1000:1 (typical)	450:1 (typical)	
	Viewing angle	±89° (typical, horizontal), ±89° (typical, vertical)	±80° (typical, horizontal), ±70° (typical, vertical)	
CD-R/CD-RW/DVD drive unit	Unit	DVD-ROM super-multi drive unit	DVD-ROM super-multi drive unit	
	Applicable disk	CD-R/CD-RW/DVD±R/DVD±R DL/DVD±RW/DVD-RAM	CD-R/CD-RW/DVD±R/DVD±R DL/DVD±RW/DVD-RAM	
	Storage capacity	8.7 GB (when DVD±R DL is used)	8.7 GB (when DVD±R DL is used)	
Hard disk drive unit	Storage capacity	500 GB (including 165 GB reserved area) Approx. 168000 images (When a 2 million-pixel image is compressed) to approx. 55000 images (When a 2 million-pixel image is not compressed)	500 GB (including 80 GB reserved area) Approx. 210000 images (When a 2 million-pixel image is compressed) to approx. 70000 images (When a 2 million-pixel image is not compressed)	
Image format		JPEG (With compression), TIFF (No compression)	JPEG/HD Photo (With compression), TIFF (No compression)	
Observable image size		20000 (H) pixels x 20000 (V) pixels (when stitched)	1600 (H) pixels x 1200 (V) pixels	
Light source	Lamp	High brightness LED	12 V, 100 W, Halogen lamp	
	Lamp life	40000 hours (reference)	1000 hours (average)	
	Colour temperature	5700K (typical)	3100K (at maximum light intensity)	
Output	Video output	DVI-I (1920 x 1080 pixels)	DVI-I (1920 x 1200 pixels)	
	Scanning frequency	Special LCD monitor	66 kHz (H), 60 Hz (V)	75 kHz (H), 60 Hz (V)
		External monitor	66 kHz (H), 60 Hz (V)	75 kHz (H), 60 Hz (V)
Input	Mouse input	USB mouse supported	USB mouse supported	
	Keyboard input	USB keyboard supported	USB keyboard supported	
	External remote input	Pause/Recording, Non-voltage input (Contact/Noncontact)	Pause/Recording, Non-voltage input (Contact/Noncontact)	
Interface	LAN	RJ-45 (10BASE-T/100BASE-TX/1000BASE-T)	RJ-45 (10BASE-T/1000BASE-T)	
	USB 2.0 Series A	6 types	8 types	
	USB 3.0 Series A	2 types	–	
Power supply	Power supply voltage	100 to 240 VAC, 50/60 Hz	100 to 240 VAC, 50/60 Hz	
	Power consumption	280 VA	340 VA	
Environmental resistance	Ambient temperature	+5 to 40°C	+5 to 40°C	
	Relative humidity	35 to 80% RH (No condensation)	35 to 80% RH (No condensation)	
Weight	Controller	Approx. 12.5 kg	Approx. 11.6 kg	
	Camera unit	Approx. 1.10 kg (VHX-5100), Approx. 1.00 kg (VHX-5020)	Approx. 0.90 kg (VHX-1020)	
	Console	Approx. 0.40 kg	Approx. 0.40 kg	
Dimensions (Excluding the projected areas)		550 (W) x 470 (H) x 200 (D) (when stored)	420 mm (W) x 416 mm (H) x 181 mm (D) (when stored)	

Basic functions: Stage

		VHX-S550E	VHX-S500E	VH-S300
XYθ stage	XY stage: Electric/Manual	Electric	Manual	Manual
	XY-motorised stage motor	2-phase stepping motor	–	–
	XY-motorised stage resolution	1 μm (typical)	–	–
	XY-motorised stage movement speed	10 mm/sec. (max.)	–	–
	XY stage moving range	±20 mm	±35 mm	±35 mm
	θ rotation angle	±90°	360°	360°
	XYθ stage size	Top surface: 171 mm x 168 mm (Centre disc: ø100)	Top surface: 190 mm x 150 mm	Top surface: 190 mm x 150 mm
Transmitted light-compatible magnification		20x or higher	–	–
Z stage	Z stage: Electric/Manual	Electric	Electric	Manual
	Z-motorised stage motor	5-phase stepping motor	5-phase stepping motor	–
	Z-motorised stage resolution	0.1 μm (typical)	0.1 μm (typical)	–
	Z-motorised stage movement speed	17 mm/sec. (max.)	17 mm/sec. (max.)	–
	Z stage moving range	49 mm	49 mm	56 mm
Ratings	Power supply voltage	100 to 240 VAC, 50/60 Hz	100 to 240 VAC, 50/60 Hz	–
	Power consumption	60 VA	60 VA	–
Environmental resistance	Ambient temperature	+5 to 40°C	+5 to 40°C	–
	Relative humidity	35 to 80% RH (No condensation)	35 to 80% RH (No condensation)	–
Weight		17.5 kg	17.0 kg	17.4 kg
Load capacity		1 kg	1 kg	1 kg

VHX-5000 (Software module details)

Software	Video recording software	Allows recording/playing back moving images.
	High quality depth composition software	Captures multiple images focused on different heights and composes a single image from them.
	Area measurement software	Measures an area of a 2D image.
	Time-lapse software	Captures images automatically at specified time intervals.
	Screen splitting software	Displays vertical, horizontal, or 4-part split screens.
	Comment input software	Allows inputting and displaying comments such as characters and markers on the observation image.
	Image improvement software	Provides image processing functions for modifying images to make observation easier.

■ Other functions

Model	VHX-5000	Console compatible	VHX-700FE	Console compatible	
Various controller functions	Auto focus function	Provided	Provided	✓	
	Image stitching	Provided	-		
	3D image stitching	Provided	-		
	High resolution image capture	Provided	-		
	Z-axis automatic stage control function	Provided	✓	Provided	✓
	One-push quick 3D function	Provided	✓	Provided	✓
	HDR+function	Provided	✓	-	
	Side album function	Provided		Provided	
	Capture condition reproduction function	Provided		Provided	
	High quality depth composition	Provided		Provided	
	Accurate D.F.D. method 3D display function	Provided (Quick method)		Provided (Quick method)	
	3D simulated illumination function	Provided		Provided	
	3D comparison function	Provided (Combination/Comparison/Difference display mode)		Provided (Combination/Comparison/Difference display mode)	
	Real-time digital zoom	1.0x to 10.0x		1.0x to 10.0x	
	Light Shift function (Height difference enhancement)	Provided (Full, partial, lateral, dark-field, bright-field, and combination illumination modes)	✓	Provided (Full, partial, and lateral illumination modes)	✓
	e-Preview mode (9 types)	Provided (Automatically lists 9 types of image modes, allowing selection of the optimal image)	✓	Provided (Automatically lists 9 types of image modes, allowing selection of the optimal image)	✓
	Glare removal function	Provided	✓	Provided	✓
	Vivid & sharp image mode	Provided		Provided	
	Supercharge shutter function	Provided	✓	Provided	✓
	Edge enhancement function	Provided (200 steps), moving images supported		Provided (200 steps), moving images supported	
	Gamma correcting function	Provided		Provided	
	Camera-shake correcting function	Provided (Moving images supported)	✓	Provided (Moving images supported)	✓
	Split function	Vertical, horizontal, 4-part, and 9-part split and combination display		Vertical, horizontal and 4-part split	
	Video recording/playback function	50 frames/sec. max. (Image size: 1600 x 1200, 800 x 600, 640 x 480)		28 frames/sec. max. (Image size: 1600 x 1200, 800 x 600, 640 x 480)	
	Timer capture function	Provided		Provided	
	Automatic unit S15 control function	Provided		Provided	
	Eucentric setting function	Provides a guide for eucentric position observation.		-	
	Real-time depth composition function	Ensures constantly focused, high depth-of-field image.	✓	-	
	High resolution HDR function	Displays a high resolution and high gradation image.	✓	-	
	High resolution observation function	Displays a high resolution image based on pixel shift technology.		-	
Simple mode	Showing a group of functions selected according to the purpose.	✓	-		
Measuring functions	TRIPLE'R function	Provided (Automatic lens connection/lens type/magnification recognition function)		Not provided. A cable is required for lens connection.	
	High-resolution dimensional measurement function	Provided		Not provided	
	Distance, angle, radius, area, and other measurement functions	Various functions provided		Various functions provided	
	Automatic count and area measurement function	Provided (Enables distance/area measurement through brightness/colour extraction)		Provided (Enables distance/area measurement through brightness/colour extraction)	
	Scale display	Various scales provided	✓	Various scales provided	✓
	Automatic edge detection	Provided		Provided	
	Auto calibration	Full-auto (Numerical input is not required)		Full-auto (Numerical input is not required)	
	One push calibration function	Provided	✓	-	
	Measurement point replacement function	Provided		Provided	
	Measurement free display function	Provided		Provided	
	Specified dimension display function	Provided		Provided	
	Measurement auxiliary function	Provided (Automatic edge extraction, multi-point input)		Provided (Automatic edge extraction, multi-point input)	
	CSV storage	Provided		Provided	
	3D height colour/scale display function	Provided (Enables X/Y/Z-axis height scale display and colour bar display related to height)		Provided (Enables X/Y/Z-axis height scale display and colour bar display related to height)	
	Height between two points measurement function	Provided		Provided	
Manual XY measurement system	XY stage measurement	Provided		Provided	
	Wide image display function	Provided		Provided	
Measuring functions (Optional functions of VHX-H4M/VHX-S15) ⁶	3D profile measurement	Provided (Displays height profile on a specified line on the 3D screen.)		Provided (Displays height profile on a specified line on the 3D screen.)	
	3D cross section profile measurement	Provided		Provided	
	3D volume measurement	Provided		Provided	
Utility	Complete style covering Observation, Recording and Measurement	All-in-one system that enables all operations for Observation, Recording, and Measurement without using a PC		All-in-one system that enables all operations for Observation, Recording and Measurement without using a PC	
	Filing system	Provided		Provided	
	Bayonet-type attachment	Provided		Provided	
	Keyboard entry	Enabled		Enabled	
	Compatible with a foot switch	Enabled		Enabled	
	User settings	Provided		Provided	
	PC mode	Provided (System protection setting available)		Provided (System protection setting available)	
Accompanying software (Free of charge, no copy restriction) (PC software)	PC communication software	Image data transfer between the VHX and PC can be performed easily. (LAN)		Image data transfer between the VHX and PC can be performed easily. (LAN)	
	3D reproduction software for the PC (Available free of charge)	The PC can reproduce a 3D image saved in the VHX.		The PC can reproduce a 3D image saved in the VHX.	
	3D HDR playback/measurement/stitched image playback software (Available free of charge)	Allows adjustment of HDR parameters and display/measurement of stitched images.		Allows measurement on the PC.	
	One-click measurement compilation software (Available free of charge)	Compiles the result of one-click measurement and transfers it to Excel [®] .		-	

1. Provides superior resolution and colour reproduction to the normal mode.
2. Provides superior colour reproduction to the high resolution mode.
3. Supported only when the multi-scan camera VHX-5100 is used.
4. Approximately 16,770,000 colours are realised with the FRC processing (dithering processing for the VHX-700FE) of the display controller.
5. The LCD monitor provided in the VHX Series is based on extremely advanced technology.
Rarely, an unlit pixel (black spot) or lit pixel (bright spot) may exist on the monitor screen. However, this is not an indication of the LCD monitor being defective.
6. The VHX-H3M is required for the VHX-700FE.

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