



IX51



## Practical Convenience and Quality Performance now enhanced with UIS2 optics



# Excellent optical performance and versatility with uncompromised quality

The IX51 continues the Olympus inverted microscope tradition by providing a very stable, compact platform ready to accept a wide range of accessories and applications. Practical convenience and quality performance are matched by excellent cost efficiency. The slim frame design enables easy attachment of multiple accessories for tissue culture and live cell imaging.



## Superb optics

#### Easy to use phase contrast system

• The 30W illumination pillar includes a condenser (N.A. 0.3, W.D. 72mm) with detachable front lens for an ultra long 182mm working distance, accommodating standing laboratory flasks or roller bottles.


#### Pre-centered phase contrast slider for quick, adjustment-free observation

• With the IX2-SLP pre-centered phase contrast slider, there is no need to center the objective each time the magnification is changed. Additionally, a common condenser annulus for 10X, 20X and 40X objectives eliminates the need to change the slider position. The result is quick and easy phase contrast observations without the need for optical adjustments and efficient throughput of routine tasks. Since the phase plate transmission rates are dictated by the objective in use, image

brightness remains constant even when the magnification is changed.


\*The IX2-SL with centering mechanism enables maximized phase contrast performance with precise optical adjustments.





#### PHC type objective minimally affected by surface tension around the container edge (CPLN10XPH, CPLFLN10XPH, LČACHN20XPH)

• The PHC type 10X and 20X objectives are designed for sharp phase contrast even at the edge of the culture dish, where the meniscus of the culture fluid can compromise other objectives. Easier multi-well observation is one of several benefits. Combined with the improved field flatness resulting from the use of UIS2 optics, this feature facilitates clear observation of cells even at the edges of a dish or well.









## Slim, compact & stable

#### Plenty of workspace

- The slim, stable, compact design with reduced footprint provides ample bench space for accessories. Shortened height of 30W halogen illumination pillar provides easy installation in a clean bench environment.
- Computer aided design with a unique aluminum alloy cast frame material combine for excellent rigidity and thermal expansion performance.





## **Ergonomics**

#### Easy, convenient operation

- Features to promote fatigue-free operation include a low stage, optional glass stage insert plate and objective colorband labels.
- Various observation tubes are provided to suit a wide range of users including a tilting binocular tube for use while sitting or standing.



## Side camera port

- Narrow frame allows a side mounted detector to access the primary image, without relay lenses, providing cameras with the highest quality image, free of aberrations.
- U-DPCAD dual camera port allows the attachment of two cameras, both with access to the primary image, with the possibility of spectral separation.



U-DPCAD double port mounted on the left side port



DP71 mounted on the left side port

## Live cell observation

#### Suitable for use with a wide range of containers

 Optical system includes objectives and condensers to accommodate different vessel types and thicknesses.

Stage adapters accommodate a variety of vessels from microscope slides to well plates.





#### Gliding stage/IX2-GS

• Designed for applications such as *C. elegans* observation, this stage moves with the touch of a finger, rotating a full 360° with 20mm X-Y travel. The flat design allows for horizontally mounted micromanipulators.

#### Plain stage (IX2-SP)+ Mechanical stage (IX-MVR)

 A mechanical stage provides X-Y movement for a variety of sample vessels. Coaxial control knobs are positioned low for user comfort.





## Cell viability

### Stable, constant environment for long-term observations

- External power supply reduces thermal expansion and electrical noise.
- All detents and click stops are designed to minimize vibration.
- For cell viability observations, CO2 incubators, chambers and a thermoplate are available.







## Fluorescence

#### **Excellent fluorescence performance**

- High throughput illuminator is 20% brighter than previous models.
- Optional L-Shaped Illuminator includes aperture and field stops and convenient access to burner centration.
- Fluorescence mirror units are designed to reduce stray light and image flare.
- Six position filter turret for maximum system flexibility.
- Improved filter performance with high efficiency coatings designed for low background and minimal crossover.









## Accessories

- An illumination pillar for a 100W halogen lamphouse allows bright, high N.A. relief contrast and DIC observations through the use of optional condensers.
- The ultra-long working distance condenser (N.A. 0.3, W.D. 73mm) with four position turret can be combined with the 100W illumination pillar for brighter 4X to 40X phase contrast and brightfield observations.







#### **UIS2** objectives

| Ob                 | jective       | N.A. | W.D.(mm) | Remarks                   |
|--------------------|---------------|------|----------|---------------------------|
| For brightfield    | PLCN4X        | 0.10 | 18.5     |                           |
|                    | PLCN10X       | 0.25 | 10.6     |                           |
|                    | PLCN20X       | 0.40 | 1.2      |                           |
|                    | LUCPLFLN40X   | 0.60 | 2.7-4    |                           |
| For phase contrast | UPLFLN4XPH    | 0.13 | 17       | PHL (for use with IX2-SL) |
|                    | CPLN10XPH     | 0.25 | 10       | PHC (for use with IX2-SL) |
|                    | PLN10XPH      | 0.25 | 10.6     | PH1 (for use with IX2-SL) |
|                    | LCACHN20XPH   | 0.40 | 3.2      | PHC (for use with IX2-SL) |
|                    | LUCPLFLN20XPH | 0.45 | 6.6-7.8  | PH1 (for use with IX2-SL) |
|                    | LCACHN40XPH   | 0.55 | 2.2      | PH2 (for use with IX2-SL) |
|                    | LUCPLFLN40XPH | 0.6  | 3.0-4.2  | PH2 (for use with IX2-SL) |
|                    | UPLFLN4XPHP*  | 0.13 | 16.4     | For use with IX2-SLP      |
|                    | CACHN10XPHP*  | 0.25 | 8.8      | For use with IX2-SLP      |
|                    | LCACHN20XPHP* | 0.40 | 3.2      | For use with IX2-SLP      |
|                    | LCACHN40XPHP* | 0.55 | 2.2      | For use with IX2-SLP      |
| J                  | CPLN10XRC     | 0.25 | 9.7      | For use with IX2-MLWCD    |
| For RC             | LCACHN20XRC   | 0.40 | 2.8      | For use with IX2-MLWCD    |
|                    | LCACHN40XRC   | 0.55 | 1.9      | For use with IX2-MLWCD    |
|                    | UPLFLN4X      | 0.13 | 17       | U,B,G                     |
| For FL             | UPLFLN10X     | 0.30 | 10       | U,B,G, w/ND filter        |
|                    | UPLFLN20X     | 0.50 | 2.1      | U,B,G                     |
|                    | LUCPLFLN20X   | 0.45 | 6.6-7.8  | U,B,G                     |
|                    | LUCPLFLN40X   | 0.6  | 2.7-4    | U,B,G                     |

<sup>\*</sup> Pre-centering objective


All UIS2 objectives and WHN eyepieces:lead-free eco-glass

#### IX51 specifications

| Microscope body                   | Revolving nosepiece                             |                         | Sextuple, simple waterproof mechanism incorporated   |  |
|-----------------------------------|---|-------------------------|--|--|
|                                   | Focus  Left side port                           |                         | 9mm stroke (from stage surface, 7mm upward and 2mm downward), coaxial coarse and fine focusing knobs (minimum fine focus graduation: 1µm, full rotation of fine focusing knobs: 100µm), upper limit stopper, torque adjustment for coarse focusing 2-step light path selection |  |
|                                   |   |                         |  |  |
| Transmitted light illuminator     | 30WHAL illumination                             | IX2-ILL30 + U-LS30-3-5  | Detachable condenser lens system (N.A. 0.3 W.D. 72mm), single filter holder (45mm, t=11mm or less), aperture iris diaphragm adjustable   |  |
|                                   | External power supply unit                      | TL-4                    | Auto voltage selector (100V / 200V)  |  |
|                                   | 100WHAL illumination                            | IX2-ILL100 + U-LH100L-3 | Pillar tilt mechanism (30° inclination angle, with vibration reducing mechanism), Condenser holder (with 50mm stroke, swing-in/out mechanism), Field iris diaphragm adjustable, 4 filter holders (45mm, 1-65mm or less)  |  |
|                                   | External power supply unit                      | TH4-100/200             | Two versions available (100V and 200V)   |  |
| Observation tube                  | Tilting binocular tube                          | U-TBI90                 | 35-85° continuous angle adjustable (eyepoint height range : 406mm-471mm), interpupillary distance adjustable between 50-76mm, diopter adjustment function, erect image, F.N. 22  |  |
|                                   | Binocular tube                                  | U-BI90CT                | Built-in focusing telescope, interpupillary distance adjustable 50-76mm, diopter adjustment function, F.N. 2   |  |
|                                   |   | U-BI90                  | Interpupillary distance adjustable 50-76mm, diopter adjustment function, F.N. 22   |  |
|                                   | Trinocular tube                                 | U-TR30H-2 + IX-ATU      | 3 step optical path selectable (observation : straight port = 100:0, 20:80, 0:100), interpupillary distance adjustable 50-76mm, diopter adjustment function, F.N. 22   |  |
| Stage                             | Cross stage with flexible right handle          | IX2-SFR                 | 50mm(X) X 50mm(Y) stroke, stage insert plate exchangeable (ø110mm)   |  |
|                                   | Cross stage with short left handle IX-SVL2      |                         | 50mm(X) X 43mm(Y) stroke, stage insert plate exchangeable (ø110mm)   |  |
|                                   | Plain stage                                     | IX2-SP                  | 232mm(X) X 240mm(Y) stage size, stage insert plate exchangeable (ø110mm)   |  |
|                                   |   | IX-MVR                  | Mechanical stage to be used with IX2-SP, 130mm(X) X 85mm(Y) stroke   |  |
|                                   | Narrow plain stage                              | IX2-KSP                 | 160mm(X) X 240mm(Y) stage size, stage insert plate exchangeable (ø110mm)   |  |
|                                   |   | CK40-MVR                | Mechanical stage to be used with IX2-KSP, 120mm(X) X 78mm(Y) stroke  |  |
|                                   | Gliding stage                                   | IX2-GS                  | Upper circular stage 360° rotatable, 20mm(X/Y) travel  |  |
| Condenser                         | Ultra long working distance                     | IX-ULWCD                | 4 positions for optical devices (for ø29mm), aperture iris diaphragm adjustable, N.A. 0.3 / W.D. 73mm  |  |
|                                   | Long working distance universal                 | IX2-LWUCD               | 5 positions for optical devices (3 positions for ø30mm and 2 position for ø38mm), aperture iris diaphragm adjustable, N.A. 0.55 / W.D. 27mm  |  |
|                                   | Long working distance Relief Contrast IX2-MLWCD |                         | 4 positions for optical devices (for ø50mm, optical Relief Contrast optical devices rotatable), aperture iris diaphragm adjustable, N.A. 0.5 / W.D. 45mm   |  |
| Eyepiece                          |   | WHN10X                  | High eyepoint, F.N. 22   |  |
|                                   |   | WHN10X-H                | High eyepoint, diopter adjustment function, F.N. 22  |  |
| Reflected light fluorescence unit | Fluorescence illuminator                        | IX2-RFAL                | L-shaped design with exchangeable F.S. and A.S. modules, two filter holder sliders (2 positions, ø32mm, t=6mm or less)   |  |
|                                   |   | IX2-RFA                 | Straight design with field iris diaphragm, filter holder slider (2 positions, ø32mm, t=6mm or less)  |  |
|                                   | Fluorescence cube turret                        | IX2-RFAC                | 6 positions in a rotating turret, built-in shutter 100W Hg lamp housing and transformer, or 75W Xe lamp housing and transformer  |  |


#### IX51 30W illumination pillar version dimensions





Weight: 20.5kg


Length marked with an asterisk (\*) may vary according to interpupillary distance.




Weight: 0.8kg Rated voltage: 55VA

#### IX51 100W illumination pillar version dimensions


(unit: mm)



Weight: 24.4kg Length marked with an asterisk (\*) may vary according to interpupillary distance.



Weight:2.2 kg Rated voltage: 200VA





- Main features of OLYMPUS Eco-products are as follows.

   Lead-free and arsenic-free Eco-glass for optics, such as lenses and prisms.

   Adoption of cardboard for packing materials without styrene foam for promoting the recycling.

  \* Some accessories are inapplicable.

Please visit our web site for further information: http://www.olympus.co.jp/en/eco-products/ Specifications are subject to change without any obligation on the part of the manufacturer.









**OLYMPUS CORPORATION obtains ISO9001/ISO14001** 



**ECO-PRODUCTS** 

OLYMPUS CORPORATION OLYMPUS LIFE AND MATERIAL SCIENCE EUROPA GMBH OLYMPUS AMERICA INC. sylvania 18034-0610, U.S.A. OLYMPUS SINGAPORE PTE LTD.

491B River Valley Road, #12-01/04 Valley Point Office Tow OLYMPUS AUSTRALIA PTY LTD.
31 Gilby Road, Mt. Waverley, VIC 3149, Melbourne, Australia.

OLYMPUS LATIN AMERICA, INC. OLYMPUS (BEIJING) SALES & SERVICE CO.,LTD.
12-13F. NCI Tower. A12 Jianguomenwai Avenue. Chaovang District. Beijing, 100022, China





This catalog is printed by environmentally-friendly waterless printing system with soy ink.